1	Intended use
2	Summary and explanation
3	Principle
4	Reagents
5	Caution
6	Storage and stability
7	Specimen collection
8	Methodology
9	Results
10	Limitations of procedure
11	Expected values
12	Performance characteristics
13	References
14	Summary of procedure

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For in-vitro diagnostic use

**Product Code: FK002._
 FS002._**

*Product manufactured by:
The Binding Site Ltd, PO Box 11712, Birmingham B14 4ZB, U.K.
www.bindingsite.co.uk*

*Telephone: +44 (0) 121 436 1000
Fax: +44 (0) 121 430 7061
e-mail: info@bindingsite.co.uk*

FDA (USA) Information
Analyte Name: Anti DNA Antibodies
Complexity Cat: High

1 INTENDED USE

This product (kit or substrate slides) is intended for use in the screening and titration of circulating anti double stranded DNA (dsDNA) autoantibodies. This autoantibody is a marker for the diagnosis and monitoring of systemic lupus erythematosus (SLE). For slides sold separately, the suitability for use with other detection systems has not been assessed but use in such assays should not necessarily be excluded.

2 SUMMARY AND EXPLANATION

Deoxyribonucleic acid is found naturally in at least three different forms: a right handed double stranded conformation (dsDNA); a left-handed double stranded conformation (Z-DNA) and a single stranded conformation (ssDNA). All forms contain determinants against which autoantibodies can be produced. dsDNA and ssDNA autoantibodies are routinely encountered in patients with autoimmune disease.

Antibodies to ssDNA (to the nucleotide bases of the DNA molecule) are not specific to a particular disease, being found in patients with rheumatoid arthritis, SLE, systemic sclerosis and other non-immunological disorders. They are a major constituent of most anti nuclear antibodies (ANAs), seen as a homogeneous pattern in indirect immunofluorescence assays.

Antibodies to dsDNA are against the phosphate deoxyribose backbone of the DNA molecule. They are less common than antibodies to ssDNA but are more clinically significant for two main reasons. Firstly, high titres of anti dsDNA are virtually only ever found in SLE patients and secondly, antibody titres correlate well with disease activity. This provides an important marker for both progression and severity of the disease and the degree of response to any subsequent treatment.

Several methods have been used over the years to detect antibodies to dsDNA, such as passive agglutination, complement fixation and RIA. The haemoflagellate, *C. luciliae* used as a substrate in indirect immunofluorescence assays, is a sensitive and specific screening test. *C. luciliae* possess an organelle called a kinetoplast; a giant mitochondrion containing dsDNA (ref. 1) but are apparently free of histones and other mammalian nuclear antigens. The advantage of *C. luciliae* based screening for dsDNA antibodies is its specificity, due to the composition of the kinetoplast (refs 2 and 3).

3 PRINCIPLE

This product utilises an indirect immunofluorescence technique (ref. 4). Patient samples and appropriate controls are incubated with the *C. luciliae* substrate. The unreacted antibodies are washed off and then an appropriate fluorescein labelled conjugate is applied. Unbound conjugate is washed off as before. Slides are viewed with a fluorescence microscope and dsDNA positive samples produce apple-green fluorescence of the kinetoplast.

4 REAGENTS

4.1 *C. luciliae* substrate slides (5- or 10-well).

Kits only

4.2 Positive control serum giving kinetoplast staining, containing 0.099% azide. Prediluted, ready for use.

4.3 Negative control serum, containing 0.099% sodium azide. Prediluted, ready for use.

4.4 Affinity purified sheep anti-human IgG (γ chain) fluorescein isothiocyanate conjugate (FITC) containing 0.099% sodium azide. Prediluted, ready for use.

- 4.5 1% Evans Blue, as an optional counterstain.
- 4.6 Phosphate buffered saline (PBS), a 20-fold concentrate in liquid form.
- 4.7 Blotters to blot around the wells after washing.
- 4.8 Mounting medium, containing an anti-fading agent (DABCO, 1, 4, diazabicyclo [2.2.2] octane).
- 4.9 Coverslips.

5 CAUTION

All donors of human serum supplied in this kit have been serum tested and found to be negative for Hepatitis B surface antigen and antibodies to Hepatitis C virus and Human Immunodeficiency Virus (HIV 1 & 2). However these tests cannot guarantee the absence of infective agents. Proper handling and disposal methods should be established as for all potentially infective material and only personnel adequately trained in such methods should be permitted to perform the procedures.

This product should only be used by suitably trained persons for the purposes stated. Adherence to the given procedure is recommended.

The kit components contain 0.099% sodium azide as a preservative and must be handled with caution – do not ingest or allow contact with skin or mucous membranes. If contact does occur wash with large volume of water and seek medical advice. Explosive metal azides may be formed with lead and copper plumbing; on disposal of reagent, flush with a large volume of water to prevent azide build up.

6 STORAGE AND STABILITY

Unopened kits/slides should be stored at 2-8°C and can be used until the expiry date given on the box label. DO NOT FREEZE. Once slides are removed from a foil bag, they should be used immediately. Diluted PBS buffer can be stored for up to one month at 2-8°C. The kit FITC conjugate should be kept out of sunlight, fluorescent or U.V. light whenever possible. All reagents should be stored at 2-8°C.

7 SPECIMEN COLLECTION

Use serum samples. Allow the blood to clot and separate the serum as soon as possible to prevent haemolysis. Sera may be stored at 2-8°C for up to 7 days prior to assay, or for prolonged storage, aliquoted and stored at -20°C or below (ref. 7). Do not freeze and thaw sera more than once. Avoid using lipaemic, haemolysed or microbially contaminated sera, as decreased titres or unclear staining patterns may occur.

8 METHODOLOGY

8.1 Materials provided

50T kit FK002.1

- 10 x *Crithidia* slide – 5-well
- 1 x 1mL *Positive Control for dsDNA* (prediluted)
- 1 x 1mL *Negative Control for dsDNA* (prediluted)
- 1 x 4.5mL *anti human IgG (γ) AFF FITC*
- 1 x 3mL *1% Evans Blue Counterstain*
- 2 x 60mL *PBS concentrate (x20)*
- 1 x 3mL *Mounting Medium*
- 20 x *Blotters*
- 10 x *Coverslips* (22 x 70mm)
- 1 x instruction leaflet

250T Kit FK002.2

- 25 x *Crithidia* slide – 10-well
- 1 x 1mL *Positive Control for dsDNA kits* (prediluted)
- 1 x 1mL *Negative Control for dsDNA kits* (prediluted)
- 1 x 12.5mL *anti human IgG (γ) AFF FITC*
- 1 x 3mL *1% Evans Blue Counterstain*
- 2 x 60mL *PBS concentrate (x20)*
- 1 x 10mL *Mounting Medium*
- 50 x *Blotters*
- 25 x *Coverslips* (22 x 70mm)
- 1 x instruction leaflet

FS002.1, FS002.2, FS002.3

- 10 x 5-well slides, 25 x 10-well slides, 100 x 10-well slides
- 1 x instruction leaflet

8.2 Additional materials required

- 8.2.1 Distilled water to dilute PBS concentrate.
- 8.2.2 Container for PBS buffer and plastic squeeze bottle for initial wash in PBS.
- 8.2.3 Micropipettes and disposable tips to apply patient samples.

- 8.2.4 Humid chamber for incubation steps (e.g. Magnetic Staining Chamber, Code BD010).
- 8.2.5 Fluorescence microscope with 495nm exciter filter and 515nm barrier filter.
- 8.2.6 If using substrate slides only, the following materials will also be required (Binding Site codes given where applicable): positive control (FA109), negative control (CON 92), anti Human IgG (γ) Aff FITC (FA004, prediluted ready for use, or AF004, dilute 1/100-1/400 for use), PBS (CON 3.3), Evans Blue (CON 93, optional), mounting medium (CON 195.1/2), blotters, coverslips.

8.3 Test procedure

Quality control

Positive and negative controls should be used every time a batch of samples is tested. It is recommended that a control is titred with every assay.

- 8.3.1 Mounting medium: Remove the mounting medium from the fridge to allow it to reach room temperature (18-28°C) before it is needed.
- 8.3.2 Dilute PBS concentrate. Dilute PBS concentrate with distilled water (1 part PBS concentrate + 19 parts distilled water) and mix. The PBS is used for diluting patient samples and as a wash buffer. N.B. Only make up total amount of kit PBS if entire kit is to be used within one month.
- 8.3.3 Dilute patient samples.
Screening: Dilute patient samples 1/10 by adding 50μL of serum to 450μL of PBS buffer.
Titration: Make serial dilutions of positive screened IgG samples with PBS buffer (e.g. 1/20, 1/40, 1/80, 1/160 and 1/320 etc.).
 For example: Take 100μL of the 1/10 dilution, mix with 100μL of PBS to give a 1/20 dilution. Repeat for further dilutions.

- 8.3.4 Substrate slides. Allow substrate slide(s) to reach room temperature (18-28°C) prior to removal from pouch(es). Label slides appropriately, place in the humid chamber and add one drop of positive and negative control to wells 1 and 2 respectively. Add 25μL of diluted patient samples to the remaining wells.

- 8.3.5 Slide incubation. Incubate slides for 30 minutes in a humid chamber at room temperature (18-28°C).

- 8.3.6 PBS wash. Remove slides from humid chamber and rinse briefly with PBS squeeze bottle. Do not squirt directly on to the wells. Place slides in a rack and immerse in PBS and agitate or stir for 10 minutes.

- 8.3.7 Addition of fluorescent conjugate. Shake off excess PBS and blot around wells using blotters provided. Return slides to humid chamber and immediately cover each well with a drop of fluorescent conjugate. DO NOT LEAVE WELLS UNCOVERED FOR LONGER THAN 15 SECONDS. Drying out of the substrate seriously affects the results.

- 8.3.8 Slide incubation. Incubate slides for 30 minutes in a humid chamber at room temperature (18-28°C) in the dark.

- 8.3.9 PBS wash. Wash again as described in step 6. Optional counterstain: Add 2-3 drops of 1% Evans Blue per 100mL of PBS prior to slide immersion.

- 8.3.10 Mounting with coverslip. Remove one slide at a time from PBS wash. Quickly dry around the wells using blotters and add a drop of mounting medium to each well. Carefully lower the slide onto the coverslip, avoiding air bubbles, but if present do not attempt to remove. Carefully wipe up any excess mounting medium.

- 8.3.11 View slides under fluorescence microscope. Finished slides should be viewed as soon as possible, but may be stored at 2-8°C in the dark for up to one week, without significant loss of fluorescence.

9 RESULTS

9.1 Quality control

The positive control should give specific apple green kinetoplast staining with or without staining of the nucleus. The negative control should show dull green staining of the organism with no discernible fluorescence. If the controls do not appear as described the assay should be repeated.


9.2 Interpretation of results

9.2.1 Negative

A sample is considered negative if specific kinetoplast staining is equivalent to or less than the negative control well, even if other structures such as the basal body, flagellum or nucleus stain.

9.2.2 Positive

A sample is considered positive if the kinetoplast staining is greater than the negative control. To clearly differentiate the kinetoplast from the nucleus, view the positive control well. The kinetoplast will always be located nearer the flagellum (see illustration below).

NB: Each laboratory should establish at which point a positive result is considered clinically significant.

INTERPRET RESULTS USING ONLY SINGLE, WELL-DEFINED ORGANISMS.

10 LIMITATIONS OF THE PROCEDURE

- 10.1 The light source, filters and optics of different makes of fluorescence microscope will influence test kit sensitivity. The performance of the microscope is significantly influenced by correct maintenance especially centring of the mercury vapour lamp and changing of the lamp after the recommended period of time.
- 10.2 Drug induced SLE patients may have dsDNA antibodies (ref. 5).
- 10.3 Patients undergoing steroid therapy may give negative results for dsDNA antibody (ref. 6).
- 10.4 Microbially contaminated serum should not be used, as false positive or negative results may occur.
- 10.5 If there appears to be an unusual amount of background staining or there is any difficulty in interpreting the results, the user may want to run a PBS control.
- 10.6 This kit is used to aid diagnosis only. A positive result suggests certain diseases which must be confirmed by clinical findings and other serological tests. The results obtained from this assay are not diagnostic proof of the presence or absence of disease.

FDA (USA) information – see front page. Slides are considered to be “Analyte Specific Reagents”. Except as a component of the dsDNA kits, analyte and performance characteristics are not established.

11 EXPECTED VALUES

C. luciliae substrate was used to test a variety of connective tissue disease patients and also 50 random normal blood donors. The results are given in the table below:

Patient group	Number	Positive <i>C. luciliae</i> staining
SLE	100	32
Rheumatoid Arthritis	10	0
Scleroderma	10	0
Sjögren's Syndrome	5	0
Normals	50	0

12 PERFORMANCE CHARACTERISTICS

12.1 Precision

Five slides from the same lot were tested with a known control sample which was diluted out to end-point. All five slides gave an end-point titre of 1/5120, demonstrating that there is no significant difference in performance between slides in the same batch. Three slides from different batches were tested with a known control sample which was diluted 1/25. All three slides gave an equivalent intensity of staining (2+).

12.2 Sensitivity

The sensitivity of the tests and the limit of detection are dependent on the microscope used. There is no international calibrator available for standardisation.

12.3 Specificity

46 sera known to contain high levels of antibodies to dsDNA and 50 normal sera were tested on The Binding Site *C. luciliae* kit and on an equivalent commercially available IFA kit. Staining was classified as positive (+), strong positive (2+) and negative (-):

Comparison kit	The Binding Site		
	2+	+	-
	23	8	0
	1	14	0
	0	0	50

13 REFERENCES

1. Aarden, L.A.; de Groot, E. R.; Felkamp, T.E.W. (1975): Immunology of DNA III. *Crithidia luciliae*, a sample substrate for the detection of anti-dsDNA with immunofluorescence techniques. Ann. NY Acad. Sci. 254: 505-515.
2. Somerfield, S.D.; Roberts, M.W.; Booth, R.J. (1981): Double-stranded DNA antibodies: comparison of four methods of detection. J. Clin. Path. 34: 1032-1035.
3. Sontheimer, R.D.; Gilliam, J.D. (1978): An immunofluorescence assay for double stranded DNA antibodies using the *Crithidia luciliae* kinetoplast as a double stranded DNA substrate. J. Lab. Clin. Med. 91: 550-558.
4. Weller, T.H.; Coons, A.H. (1954): Fluorescent antibody studies with agent of Varicella and Herpes Zoster propagated in Vitro: Proc. Soc. Exp. Biol. Med. 86: 789-794.
5. Epstein, W.V. (1975): Specificity of SLE serum antibody for single-stranded and double-stranded DNA configuration. J. Rheum. 2: 215-220.
6. Ballou, S.P.; Kushner, I. (1979): Anti-native DNA detection by the *Crithidia luciliae* method. Arthritis Rheum. 22: 321-328.
7. Protein Reference Handbook of Autoimmunity (3rd Edition) 2004 Ed A Milford Ward, J Sheldon, GD Wild. Publ. PRU Publications, Sheffield. 14.

14 SUMMARY OF PROCEDURE

- 14.1 Equilibrate mounting medium to room temperature.
- 14.2 Dilute PBS with distilled water.
- 14.3 Dilute patient sera 1/10 with PBS.
- 14.4 Remove slides from refrigerator and equilibrate to room temperature (18-28°C).
- 14.5 Remove slide from foil bag and place in a humid chamber. Add 25µL of controls and patient sera.
- 14.6 Incubate for 30 minutes at room temperature (18-28°C).
- 14.7 Rinse slides with a stream of PBS. Wash slides for 10 minutes in a rack.
- 14.8 Blot around each well, return slide to humid chamber and immediately add a drop of fluorescent conjugate.
- 14.9 Incubate slides for 30 minutes.
- 14.10 Wash again as described in step 7.
- 14.11 Blot around wells; add mounting media to each well and coverslip.
- 14.12 View slides under a fluorescence microscope.

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13 REFERENZEN

- L. A. Aarden; E. R. deGroot; T. E. W. Felkamp (1975): "Immunology of DNA III. *Crithidia luciliae* a sample substrate for the detection of anti-dsDNA with immunofluorescence technique." Ann. NY Acad. Sci. 254, 505 – 515.
- S. D. Sommerfield; M. W. Roberts; R. J. Booth (1981): "Double-stranded DNA antibodies: comparison of four methods of detection." J. Clin. Path. 34, 1032 – 1035.
- R. D. Sontheimer; J. D. Gilliam (1981): "An immunofluorescence assay for double stranded DNA antibodies using the *Crithidia luciliae* kinetoplast as a double stranded DNA substrate." J. Lab. Clin. Med. 91, 550 – 558.
- T. H. Weller; A. H. Coons (1954): "Fluorescent antibody studies with agent of Varicella and Herpes Zoster propagated in Vitro." Proc. Soc. Exp. Biol. Med. 86, 789 – 794.
- W. V. Epstein (1975): "Specificity of SLE serum antibody for single-stranded and double-stranded DNA configuration." J. Rheum. 2, 215 – 220.
- S. P. Ballou; I. Kushner (1979): "Anti-native DNA detection by the *Crithidia luciliae* method." Arthritis Rheum. 22, 321 – 328.
- Protein Reference Handbook of Autoimmunity (3rd Edition) 2004 Ed A Milford Ward, J Sheldon, GD Wild. Publ. PRU Publications, Sheffield. 14.

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13 BIBLIOGRAPHIE

1. Aarden, L.A.; de Groot, E. R.; Felkamp, T.E.W. (1975): Immunology of DNA III. *Crithidia luciliae*, a simple substrate for the detection of anti-dsDNA with immunofluorescence techniques. Ann. NY Acad. Sci. 254: 505-515.
2. Somerfield, S.D.; Roberts, M.W.; Booth, R.J. (1981): Double-stranded DNA antibodies : comparison of four methods of detection. J. Clin. Path. 34: 1032-1035.
3. Sontheimer, R.D.; Gilliam, J.D. (1978): An immunofluorescence assay for double stranded DNA antibodies using the *Crithidia luciliae* kinetoplast as a double stranded DNA substrate. J. Lab. Clin. Med. 91: 550-558.
4. Weller, T.H.; Coons, A.H (1954): Fluorescent antibody studies with agent of Varicella and Herpes Zoster propagated in vitro : Proc. Soc. Exp. Biol. Med. 86: 789-794.
5. Epstein, W.V. (1975): Specificity of SLE serum antibody for single-stranded and double-stranded DNA configuration. J. Rheum. 2: 215-220, 1975.
6. Ballou, S.P.; Kushner, I. (1979): Anti-native DNA detection by the *Crithidia luciliae* method. Arthritis Rheum. 22: 321-328.
7. Protein Reference Handbook of Autoimmunity (3rd Edition) 2004 Ed A Milford Ward, J Sheldon, GD Wild. Publ. PRU Publications, Sheffield. 14.

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1. Aarden, L.A.; de Groot, E. R.; Felkamp, T.E.W. (1975): Immunology of DNA III. *Crithidia luciliae*, a sample substrate for the detection of anti-dsDNA with immunofluorescence techniques. Ann. NY Acad. Sci. 254: 505-515.
2. Somerfield, S.D.; Roberts, M.W.; Booth, R.J. (1981): Double-stranded DNA antibodies: comparison of four methods of detection. J. Clin. Path. 34: 1032-1035.
3. Sontheimer, R.D.; Gilliam, J.D. (1978): An immunofluorescence assay for double stranded DNA antibodies using the *Crithidia luciliae* kinetoplast as a double stranded DNA substrate. J. Lab. Clin. Med. 91: 550-558.
4. Weller, T.H.; Coons, A.H. (1954): Fluorescent antibody studies with agent of Varicella and Herpes Zoster propagated in Vitro: Proc. Soc. Agent. Biol. Med. 86: 789-794.
5. Epstein, W.V. (1975): Specificity of SLE serum antibody for single-stranded and double-stranded DNA configuration. J. Rheum. 2: 215-220.
6. Ballou, S.P.; Kushner, I. (1979): Anti-native DNA detection by the *Crithidia luciliae* method. Arthritis Rheum. 22: 321-328.
7. Protein Reference Handbook of Autoimmunity (3rd Edition) 2004 Ed A Milford Ward, J Sheldon, GD Wild. Publ. PRU Publications, Sheffield. 14.

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13 REFERENCIAS

1. Aarden, L.A.; de Groot, E. R.; Felkamp, T.E.W. (1975): Immunology of DNA III. *Crithidia luciliae*, a sample substrate for the detection of anti-dsDNA with immunofluorescence techniques. Ann. NY Acad. Sci. 254: 505-515.
2. Somerfield, S.D.; Roberts, M.W.; Booth, R.J. (1981): Double-stranded DNA antibodies: comparison of four methods of detection. J. Clin. Path. 34: 1032-1035.
3. Sontheimer, R.D.; Gilliam, J.D. (1978): An immunofluorescence assay for double stranded DNA antibodies using the *Crithidia luciliae* kinetoplast as a double stranded DNA substrate. J. Lab. Clin. Med. 91: 550-558.
4. Weller, T.H.; Coons, A.H. (1954): Fluorescent antibody studies with agent of Varicella and Herpes Zoster propagated in Vitro: Proc. Soc. Exp. Biol. Med. 86: 789-794.
5. Epstein, W.V. (1975): Specificity of SLE serum antibody for single-stranded and double-stranded DNA configuration. J. Rheum. 2: 215-220.
6. Ballou, S.P.; Kushner, I. (1979): Anti-native DNA detection by the *Crithidia luciliae* method. Arthritis Rheum. 22: 321-328.
7. Protein Reference Handbook of Autoimmunity (3rd Edition) 2004 Ed A Milford Ward.J.Sheldon, GD Wild. Publ. PRU Publications, Sheffield. 14.

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