MONKEY STOMACH IFA SLIDES

For *in-vitro* diagnostic use

**Product Code: FS230.A
FS230.1**

Product manufactured by:
The Binding Site Ltd, PO Box 11712, Birmingham B14 4ZB, UK
www.bindingsite.co.uk

Telephone: +44 (0)121 436 1000
Fax: +44 (0)121 430 7061
e-mail: info@bindingsite.co.uk

1 INTENDED USE

This product is intended for use in the screening and titration of circulating autoantibodies in human serum as an aid in the diagnosis and treatment of various autoimmune diseases. The two major autoantibodies detected are antigastric parietal cell (AGPCA) and antismooth muscle (ASMA).

2 SUMMARY AND EXPLANATION

Although the most commonly used substrate for detection of autoantibodies is rodent tissue, it can be argued that this is not ideal for diagnostic purposes. Some clinically significant autoantibodies bind poorly to rodent antigens and, conversely, heterophile antibodies bind to rodent tissues but are not clinically relevant (refs 2 and 3).

Monkey tissues are more similar to human tissues than rodent tissues with regard to both antigenic expression and organ structure. Unfortunately, monkey and human immunoglobulins are also very similar, sharing approximately 98% of their amino acid sequences. The background fluorescence, produced by staining of endogenous immunoglobulin, is therefore similar using either human or monkey tissues with conventional anti-human immunoglobulin conjugates. For this reason, The Binding Site's monkey adsorbed anti-human immunoglobulin FITC conjugates (FA003.M, FA004.M) are recommended for use with these slides. The absence of non-specific tissue staining, with these conjugates, enhances the sensitivity and facilitates interpretation.

Particular autoantibodies are associated with a number of different diseases. SMA are frequently associated with chronic active hepatitis and primary biliary cirrhosis, but are also detectable in low concentrations in various other conditions. GPC occur in the serum of most patients with pernicious anaemia

3 PRINCIPLE

These slides are used in an indirect immunofluorescence technique where patient samples and appropriate controls are incubated with the sections. The unreacted antibodies are washed off and then appropriate fluorescein- labelled conjugates are applied. Unbound conjugate is washed off as before. Slides are viewed with a fluorescence microscope and positive samples produce apple-green fluorescence that corresponds to areas of the section where autoantibody has bound (ref. 1).

4 REAGENTS

Monkey stomach sections on 5-well slides individually wrapped in a foil pouch containing a desiccant.

5 CAUTION

Proper handling and disposal methods should be established for all potentially infective samples tested with this product; only personnel adequately trained in such methods should be permitted to perform the procedures.

6 STORAGE AND STABILITY

Unopened slides should be stored at 2-8°C and can be used until the given expiry date. DO NOT FREEZE. Once slides are removed from a foil bag, they should be used immediately.

7 SPECIMEN COLLECTION

Blood samples should be collected by venepuncture, allowed to clot naturally and the serum separated as soon as possible to prevent haemolysis. The serum may be stored at 2-8°C for up to 7 days prior to assay (ref. 4), or for prolonged storage, aliquoted and stored at -20°C or below. DO NOT freeze and thaw sera more than once. Avoid using lipaemic, haemolysed or microbially contaminated sera as decreased titres or unclear staining patterns may occur.

8 METHODOLOGY

8.1 Materials provided

FS230.A, FS230.1

- 8.1.1 1 x *Monkey stomach*– 5 well slide or 10 x *Monkey Stomach 5 well slides* respectively.
- 8.1.2 1 x instruction leaflet.

8.2 Additional materials required

- 8.2.1 PBS for sample diluent and washes.
- 8.2.2 Container for PBS buffer.
- 8.2.3 Micropipettes and disposable tips to apply patient samples.
- 8.2.4 Humid chamber for incubation steps (e.g. magnetic staining chamber, Code BD010).
- 8.2.5 Fluorescence microscope with 495nm exciter filter and 515nm barrier filter.
- 8.2.6 Plastic squeeze bottle for initial wash in PBS.
- 8.2.7 Additional components may be obtained from The Binding Site: PBS (CON 3.3), negative control (CON 92), SMA positive control (FA134), GPC positive control (FA125), anti-Human IgG (H+L) monkey adsorbed conjugate (FA003.M), 1% Evan's Blue (CON 93) and Mounting medium (CON195).

8.3 Test Procedure

Quality Control

Positive and negative controls should be used every time samples are tested.

- 8.3.1 Mounting Medium: Remove the mounting medium from the fridge to allow it to reach room temperature (18°C-28°C) before it is needed.

8.3.2 Dilute Patient Samples

Screening: Dilute patient samples 1/20 by adding 10µL of serum to 190µL of PB buffer.

Titration: Make serial dilutions of positive screened samples with PBS buffer (e.g. 1/20, 1/40, 1/80 and 1/160 etc). For example: Take 100µL of the 1/20 dilution, mix with 100µL PBS to give a 1/40 dilution (repeat for further dilutions).

- 8.3.3 Substrate Slides. Allow substrate slide(s) to reach room temperature (18-28°C) prior to removal from pouch(es). Label slides appropriately, place in the humid chamber and add positive and negative controls to appropriate wells. Add 50µL of diluted patient samples to the remaining wells.

- 8.3.4 Slide Incubation. Incubate slides for 30 minutes in a humid chamber at room temperature (18-28°C).

- 8.3.5 PBS Wash. Remove slides from humid chamber and rinse briefly with PBS squeeze bottle. Do not squirt directly on to the wells. Place slides in a rack and immerse in PBS and agitate or stir for 5-10 minutes.

- 8.3.6 Addition of fluorescent conjugate. Shake off excess PBS and blot around wells. Return slides to humid chamber and immediately cover each well with a drop of appropriately diluted fluorescent conjugate. DO NOT LEAVE WELLS UNCOVERED FOR LONGER THAN 15 SECONDS. Drying out of the substrate seriously affects the results. The use of a monkey adsorbed conjugate will greatly enhance results (e.g. FA003.M).

- 8.3.7 Slide Incubation. Incubate slides for 30 minutes in humid chamber at room temperature (18-28°C), in the dark.

- 8.3.8 PBS Wash. Wash again as described in step 8.3.5. OPTIONAL COUNTERSTAIN. Add 2-3 drops of 1% Evans Blue per 100mL of PBS prior to slide immersion.

- 8.3.9 Mounting with coverslip. Remove one slide at a time from PBS wash. Quickly dry around the wells and add a drop of mounting medium to each well. Carefully lower the slide onto the coverslip, avoiding air

bubbles but, if present, do not attempt to remove. Wipe excess medium from around edge of coverslip.

- 8.3.10 View slides under fluorescence microscope. Slides may be stored for up to 3 days at 2-8°C, in the dark, without significant loss of fluorescence.

9 RESULTS

9.1 Quality Control

- 9.1.1 A serum sample containing SMA should give bright apple green fluorescent staining of the muscularis layers and interglandular fibres of the stomach.
- 9.1.2 A serum sample containing GPC should give bright apple green fluorescent staining of the parietal cells of the stomach.
- 9.1.3 A negative control should show dull green staining in all the tissue, with no discernible fluorescence.

If controls do not appear as described, the test is invalid and should be repeated.

9.2 Interpretation of Results

See references 2 and 3 for colour photographic examples of these patterns. Results are reported as positive or negative.

N.B: Each laboratory should establish at which point a positive result is considered clinically significant.

10 LIMITATIONS OF PROCEDURE

- 10.1 The light source, filters and optics of different makes of fluorescence microscopes will influence the sensitivity of the kit. The performance of the microscope is significantly influenced by correct maintenance especially centring of the mercury vapour lamp and changing of the lamp after the recommended period of time.
- 10.2 Suitability for use with other manufacturers' IFA reagents has not been assessed but use with such reagents should not necessarily be excluded.
- 10.3 This test alone should not be considered diagnostic. All other factors including the clinical history of the patient and other serological or biopsy results must also be taken into account.

11 REFERENCES

1. Weller T. H. & Coons A. H. (1954). Fluorescent antibody studies with agent of Varicella and Herpes Zoster propagated in vitro: Proc. Soc. Exp. Biol. Med. **86**: 789-794.
2. Bradwell A. R. *et al* (1997). Atlas of autoantibody patterns on tissues. Publ. The Binding Site Ltd., Birmingham UK.
3. Bradwell A. R. *et al* (1999). Advanced atlas of autoantibody patterns. Publ. The Binding Site Ltd., Birmingham UK.
4. Protein Reference Handbook of Autoimmunity (3rd Edition) 2004. Ed. A Milford Ward, GD Wild. Publ. PRU Publications, Sheffield. 14.

12 SUMMARY OF PROCEDURE

- 12.1 Equilibrate mounting medium to room temperature.
- 12.2 Dilute PBS with distilled water.
- 12.3 Dilute patient sera 1/20 with PBS.
- 12.4 Equilibrate substrate slides to room temperature (18-28°C).
- 12.5 Apply 50-100µL positive and negative controls and diluted patient sera to appropriate wells.
- 12.6 Incubate in a humid chamber for 30 minutes.
- 12.7 Wash for 5-10 minutes in PBS.
- 12.8 Blot around each well and immediately cover each well with a drop of conjugate.
- 12.9 Incubate as in step 12.6.
- 12.10 Wash as step 12.7.
- 12.11 Mount.
- 12.12 View slide under fluorescence microscope.

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