

## PERFORMANCE OF THREE NEW CHEMILUMINESCENT ASSAYS FOR CELIAC DISEASE IN PEDIATRIC PATIENTS SHOWING DIFFICULT CLINICAL AND SEROLOGIC PRESENTATION

Marta Molero (1) BS; Carme Farré (1) PhD; Silvia Casas (2) PhD; Priscilla Carrion (2) BS; Alison Lewis (2) BS; Kitty Li (2) BS; Peter Martis (2) BS and Rufus Burlingame (2) PhD.  
 (1) Hospital Sant Joan de Déu. Barcelona. Spain; (2) INOVA Diagnostics Inc. San Diego. CA

### Introduction

The tissue transglutaminase antibody (tTG-IgA) immunoassay is the universally recommended screening test for CD. In IgA deficiency cases, testing with tTG-IgG is recommended. Because of the inferior accuracy of the antigliadin assay, this test is refused. The use of deamidated gliadin peptides (DGP) to improve specificity awaits confirmation by large-scale validation. The clinical of CD is very unspecific, from asymptomatic to classical symptoms. Small bowel biopsy is the confirmatory cornerstone test. Clinical suspicion and/or positive serological test are necessary to recommend intestinal study. The histological (serological and/or clinical) response to gluten free diet (GFD) confirms CD diagnosis.

### Objective

Determine if the new serological tests can aid in the diagnosis of CD in paediatric patients with difficult clinical and serologic presentation.

### Results

1 - Specificity: all 29 supposed non-CD paediatric patients were negative for both IgA anti-tTG and IgG anti-tTG antibodies, while one patient was positive on the DGP-screen on the BIO-FLASH™.

2 - Sensitivity: all 27 asymptomatic CD patients were positive on the IgA anti-tTG, 12/27 on the IgG anti-tTG and 25/27 on the DGP-screen.

3 - Border zone: 10 of 27 CD patients on GFD were low positive on IgA anti-tTG, all were negative on the IgG anti-tTG and 5/27 were low positive on the DGP-screen.

4 - All 3 clinical CD patients diagnosed in spite of their negative serological markers were also negative in the new IgA anti-tTG and IgG anti-tTG assays. One of them (1/3) was positive on the DGP-screen.

We highlight the extended linearity of IgA anti-tTG on the chemiluminescent BIO-FLASH™ assay that allows quantitation of a wide concentrations range.

### Material and methods

Samples were chosen from four groups of patients:

1 – Specificity: 29 supposed non-CD paediatric patients from 1 to 2 years old with digestive symptoms but negative serological markers.

2 – Sensitivity: 27 asymptomatic CD patients (Marsh III) detected by positive serological markers obtained in an at risk population.


3 – Border zone: 27 CD patients on GFD with low levels of serological markers.

4 – And finally, 3 CD patients (Marsh III) displaying classical CD symptoms in spite of negative serological markers.


All serum samples were tested on 3 new automated chemiluminescent immuno-assays in the BIO-FLASH™ Analyzer for: IgA anti-tTG antibodies, IgG anti-tTG antibodies and IgA/IgG anti-DGP antibodies (DGP-screen).

(Samples were previously analyzed by the established IgA anti-tTG EIA ImmunoCAP250 from PHADIA)

Figure 1 and 2 illustrate the good correlation between the new IgA anti-tTG run on the BIO-FLASH™ and the established IgA anti-tTG EIA run on the ImmunoCAP250 PHADIA in all patients (Figure 1) and only in border zone (Figure 2). Positive value of BIO-FLASH™ ≥ 20 U/ml. Positive value of PHADIA ≥ 3.8 U/ml.



**Figure 1. Correlation IgA anti-tTG BIO-FLASH™ and PHADIA.** These results include non-CD patients, asymptomatic CD patients and CD patients on GFD.



**Figure 2. Correlation of IgA anti-tTG BIO-FLASH™ and PHADIA in border zone.** The border zone includes 27 CD patients on a gluten free diet.

### Conclusions

- Both the new IgA anti-tTG and DGP-screen chemiluminescent immuno-assay methods run on the automated BIO-FLASH™ Analyzer are useful to detect antibodies in CD patients.
- The IgA anti-tTG on the BIO-FLASH™ is equal to the results from the IgA anti-tTG on the EIA ImmunoCAP250 PHADIA.
- The IgA anti-tTG on the BIO-FLASH™ is more efficient than the DGP-screen on the BIO-FLASH™ to detect asymptomatic CD patients in at risk populations because it detects 27 of 27 patients, while the DGP-screen finds 25 of 27 patients.
- In CD patients on GFD treatment, the DGP-screen yields less positive results than IgA anti-tTG, suggesting that the anti-DGP antibodies disappeared more quickly than IgA anti-tTG in GFD patients of our series.
- In our criteria, IgG anti-tTG analysis is only recommended in IgA deficient patients.
- The DGP-screen result was positive in 1 of 3 diagnosed CD patients who were negative for the conventional IgA anti-tTG autoantibodies